

Every OneTouch® meter is designed to make it easy to test your blood glucose and help you manage your diabetes.

Step
1

Insert a OneTouch Verio® test strip



Step
2

Apply blood to either side of the strip



Virtually pain-free testing

And that's it.

Your results appear on screen. Our ColourSure® technology will instantly show if your results are **in or out of range**.



Connect to the **OneTouch Reveal® App** for even more insight. See reverse side.

Current numerical result and healthcare professional recommendations should be considered in treatment decisions.

Always read and follow the label. © 2020-2022 LifeScan IP Holdings, LLC. Trademarks are used under license by LifeScan IP Holdings, LLC. All rights reserved. CA-VRF-2000072

ONETOUCH®

Sync data, seamlessly, from your OneTouch® meter directly to your smartphone. Let's get started....

1. Download the FREE OneTouch Reveal® app

- Visit the app store on your mobile device, search "OneTouch Reveal", then install the app



2. Activate Bluetooth® on your OneTouch® meter

For the OneTouch Verio Reflect® meter

- Turn your meter on. Press and hold **OK** until the startup screen appears
- Go to the settings menu > select and activate Bluetooth®



Bluetooth®
Don't forget to also activate Bluetooth® on your mobile device.

For the OneTouch Verio Flex® meter

- Press the **OK** button to turn your meter on
- Press the **OK** and **^** buttons at the same time to activate Bluetooth®

3. Connect

- Open the OneTouch Reveal® app on your mobile device and follow the on-screen instructions
- When prompted, enter the PIN shown on your meter screen into your mobile device and tap **Pair**



You can also view your account on your computer at OneTouchReveal.ca

12 3456

Example of a PIN number

Now you're ready to test and sync!

(Optional) Share data with your Care Team

- After setup, tap **MORE** and select **Connections**
- Tap **Add a Clinic** and enter the **Clinic Code** provided by your Care Team and follow the rest of the on-screen instructions

Enter Clinic Code

--	--	--	--	--	--

Discover more at OneTouch.ca or call **1-800-663-5521**.

Apple, App Store, and the Apple logo are trademarks of Apple Inc., registered in the US and other countries. Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by LifeScan Scotland Ltd. and its affiliates is under license. Other trademarks and trade names are those of their respective owners. Always read and follow the label. © 2020-2022 LifeScan IP Holdings, LLC. Trademarks are used under license by LifeScan IP Holdings, LLC. All rights reserved. CA-VRF-2000072

ONETOUCH®