
ONETOUCH[®] Verio[®]

Test Strips
Bandelettes de test

AW 067069008A

For *in vitro* diagnostic use. For self-testing.

IMPORTANT: Please read this insert and your User Guide or Owner's Booklet **that came with your system** before using. The User Guide or Owner's Booklet has full instructions for blood glucose and control solution testing. **Do Not** use new test strips if the vial is open or damaged in any way. This could lead to error messages or inaccurate blood glucose results. Contact Customer Service right away if you have questions about your test strips, test results or these instructions. Contact the OneTouch[®] Customer Care Line at 1 800 663-5521, Monday-Friday, 9am-8pm Eastern, 6am-5pm Pacific, or visit www.OneTouch.ca.

Intended Use

The OneTouch Verio[®] Test Strips are intended for use with the OneTouch Verio[®] Family of Blood Glucose Monitoring Systems for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertip. The OneTouch Verio[®] Test Strips and associated meters are intended for use by people with diabetes at home as an aid to monitor the effectiveness of diabetes control.

The OneTouch Verio[®] Family of Blood Glucose Monitoring Systems are intended to be used by a single patient and should not be shared.

Storage and Handling

- Store the test strip vial in a cool dry place between 5°C and 30°C. **Do Not** refrigerate. Keep away from direct sunlight and heat. Exposure to temperatures and/or humidity outside the storage limits may cause inaccurate readings.
- The expiration date is printed on the vial.
- The discard date is 6 months after you first open a test strip vial. When you open a new vial, write the discard date on the label.
- **Do Not** use after the expiration or discard date, whichever comes first.
- Store your test strips in the **vial only**. To avoid damage or contamination, **Do Not** put test strips in any other container.
- Open the **vial only** when taking out a test strip for use.
- Close the vial lid tightly right away after removing a test strip. Use each test strip **as soon as you take it out** of the vial.
- **Do Not** use test strips from a vial that is damaged or left open to air.
- **Do Not** get dirt, food or liquids on the test strip. You may touch the test strip anywhere on its surface with clean, dry hands.
- **Do Not** bend, cut or alter test strips in any way.
- Test strips are for single use only.
- **Never reuse a test strip that already has blood or control solution on it.**
- Make sure your meter and test strips are about the same temperature before you test.
- Apply only OneTouch Verio[®] Control Solution or a blood sample to the test strip.
- **Do Not** put the used test strip back in the vial after taking a test.
- Used test strips may be considered biohazardous waste in your area. Be sure to follow your healthcare professional's advice or local laws when throwing them away.

⚠️WARNING: Keep the test strip vial away from children. Test strips are a choking hazard. **Do Not swallow test strips. The test strip vial may contain drying agents that are harmful if inhaled or swallowed and may cause skin or eye irritation. Do Not ingest or swallow any items.**

Testing your blood glucose

Read your system User Guide or Owner's Booklet for instructions on running a blood test and getting a blood sample.

Test Results

The OneTouch Verio[®] Family of Meters shows results between 1.1 and 33.3 mmol/L.

Read these cautions when your test results are higher or lower than what you expect.

⚠️CAUTION:

Low glucose results

If your blood glucose result is below 4.0 mmol/L or a low glucose warning is displayed, (meaning the result is less than 1.1 mmol/L), it may mean hypoglycemia (low blood glucose). You may need to treat low glucose right away. Follow your healthcare professional's advice. The result could be due to a test error, but it is safer to treat first, then test again.

Dehydration and low glucose results

You may get false low glucose results if you are severely dehydrated. If you think this is the case, contact your healthcare professional right away.

High glucose results

If your test result is above 10.0 mmol/L, it may mean high blood glucose (hyperglycemia). You should think about re-testing. Talk to your healthcare professional if you are concerned about high blood glucose.

A high glucose warning is shown when your blood glucose level is over 33.3 mmol/L. You may have very high blood glucose (severe hyperglycemia). Re-test your blood glucose. If the result is a high glucose warning again, this indicates a severe problem with your blood glucose control. Obtain and follow instructions from your healthcare professional right away.

Repeated unexpected glucose results

If you get repeated unexpected results:

- Run a control solution test to check your system.
- Make sure you follow all instructions in your User Guide or Owner's Booklet.
- If your symptoms don't match your blood glucose results, call your healthcare professional.

Never ignore symptoms or make significant changes to your diabetes management program without speaking to your healthcare professional.

Range of Expected Values

Blood glucose management requires the help of a healthcare professional. Talk to your healthcare professional about:

- The target blood glucose range that is right for you.
- When to test.
- What your blood glucose results mean.

Expected blood glucose levels for non-pregnant people without diabetes¹:

Time	Range, mmol/L
Fasting	Less than 7.0
2 hours after meals	Less than 11.1

¹Diabetes Canada Clinical Practice Guidelines Expert Committee. *Diabetes Canada 2018 Clinical Practice Guidelines for the Prevention and Management of Diabetes in Canada*. Can J Diabetes. 2018;42(Suppl 1):S1-S325.

Checking the System

Use only OneTouch Verio[®] Control Solution

A control solution test is used to check that:

- The meter and test strips are working together properly.
- You are doing the test correctly.

Do a control solution test

- When you open a new vial of test strips.
- If you think the meter or test strips are not working properly.
- If you get repeated unexpected blood glucose results.
- If you drop or damage the meter.

Read your system User Guide or Owner's Booklet for instructions on how to run a control solution test.

Limitations of Procedure

OneTouch Verio[®] Test Strips give accurate results under the following conditions:

- Use only fresh capillary whole blood. **Do Not** use serum or plasma.
- Hematocrit is the percentage of red blood cells in the blood. Extremes in hematocrit may affect test results. Hematocrit levels below 20% may cause false high readings. Hematocrit levels over 60% may cause false low readings. Talk to your healthcare professional if you don't know your hematocrit level.
- OneTouch Verio[®] Test Strips may be used at altitudes up to 3048 meters (10,000 feet) with no effect on test results.
- Interferences: Uric acid and other reducing substances (when occurring in blood at high normal or at high therapeutic concentrations) do not significantly affect results. However, abnormally high concentrations in blood may cause inaccurately low results.
- Test results may be falsely low if the patient is severely dehydrated, in shock or in a hyperosmolar state (with or without ketosis).
- Critically ill patients should not be tested with blood glucose meters.

⚠️CAUTION:

- The OneTouch Verio[®] Test Strips should not be used for patients within 24 hours of receiving a D-xylose absorption test as it may cause inaccurately high results.
- **Do Not** use the OneTouch Verio[®] Family of Blood Glucose Monitoring Systems when PAM (Pralidoxime) is known or suspected to be in the patient's whole blood sample.

Test Principle

The OneTouch Verio[®] Family of Meters is plasma-calibrated for easy comparison of results with lab methods. Glucose in the blood sample mixes with special chemicals on the test strip to produce an electrical current. The strength of the current changes with the amount of glucose in the blood sample. The OneTouch Verio[®] Family of Meters measures the current, calculates your blood glucose level and shows your result.

Reagent Composition

Each test strip contains: flavin adenine dinucleotide dependent glucose dehydrogenase or FAD-GDH (from *Aspergillus* sp.) - 2 U; potassium ferricyanide 41 µg; and other ingredients (buffer, etc.). The vial contains a drying agent.

Performance Characteristics

See your system User Guide or Owner's Booklet.

IMPORTANT: For instructions, technical data and a list of symbols used, read your system User Guide or Owner's Booklet.

OUR COMMITMENT TO YOU

Our goal is to provide you with quality healthcare products and dedicated customer service. If you are not fully satisfied with this product, or if you have questions about the use of any LifeScan product, contact Customer Service. Contact the OneTouch[®] Customer Care Line at 1 800 663-5521, Monday-Friday, 9am-8pm Eastern, 6am-5pm Pacific, or visit www.OneTouch.ca. If you cannot reach Customer Service, contact your healthcare professional for advice.

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Description des symboles

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Contents covered by one or more of the following U.S. patents: 6,179,979, 6,716,577, 6,749,887, 7,045,046 and 7,846,312. Use of these test strips and associated monitoring device is protected under the following U.S. patents: 6,193,873, 6,475,372, 6,797,150, 7,498,132, and 8,449,740. Purchase of the associated monitoring device does not act to grant a use license under these patents. Such a license is granted only when the associated monitoring device is used with OneTouch Verio® Test Strip. No test strip supplier other than LifeScan, Inc. is authorized to grant such a license. The accuracy of results generated with LifeScan meters using test strips manufactured by anyone other than LifeScan has not been evaluated by LifeScan.

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Read instructions.
Lisez les instructions.

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© 2014–2018 LifeScan, Inc.
LifeScan Europe GmbH