

## How to test guide

The OneTouch Verio Flex™ meter takes the guesswork out of your numbers...

ColourSure™  
TECHNOLOGY



**Step 1:** Insert a OneTouch Verio® test strip

+



**Step 2:** Apply blood to either side of the strip

=



**And that's it.** Your results appear on screen with an arrow that shows if your results are **Low**, **In Range**, or **High**

Before getting your blood glucose results, please read the OneTouch Verio Flex™ meter owner's booklet carefully.

...with simple, 2-step testing

All trademarks are used under licence by LifeScan Canada.  
© 2016 LifeScan Canada. All rights reserved. NACO/VFX/1116/0517

**ONETOUCH**®  
every touch is a step forward™

## Comment effectuer un test

Le lecteur OneTouch Verio Flex<sup>MC</sup> élimine l'incertitude concernant vos résultats de glycémie...

TECHNOLOGIE  
ColourSure<sup>MC</sup>



**1<sup>re</sup> étape :** Insérez une bandelette de test OneTouch Verio<sup>®</sup>

+



**2<sup>e</sup> étape :** Appliquez le sang de n'importe quel côté de la bandelette

=



**Le tour est joué.** Vos résultats s'affichent à l'écran avec une flèche qui indique si votre glycémie est **inférieure à l'objectif**, **dans l'objectif** ou **supérieure à l'objectif**.

Avant d'effectuer un test de glycémie, veuillez lire attentivement le Manuel d'utilisation du lecteur OneTouch Verio Flex<sup>MC</sup>.

...car il permet d'effectuer un test en deux étapes simples.